coltène whaledent®

ParaPost® TAPER LUX™
Esthetic Post System

Instructions for Use

CE
0344

EN **Instructions For Use** 3
DE **Gebrauchsanleitung** 7
FR **Mode d'emploi** 11
ES **Instrucciones de uso** 15
IT **Istruzioni per l'uso** 19
NL **Gebruiksaanwijzing** 23
SV **Användarinstruktioner** 27
DA **Brugsanvisning** 31
FI **Käyttöohjeet** 35
NO **Bruksanvisning** 38
PT **Instruções para uso** 41
EL **Οδηγίες Χρήσης** 45

One-Office-Visit Technique

For illustrations see page 50

Caution: Federal law restricts this device to sale by or on the order of a dentist.

We urge the use of Hygenic® Rubber Dam with this procedure.

1. After endodontic therapy has been completed, prepare the tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect **(Fig. 01)**. *If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.*
2. Use a radiograph to determine the appropriate diameter and depth of the planned post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)* to prevent perforation or weakening of root walls, and at least 4 - 5 mm of gutta percha must be left intact to protect the apical seal.
3. Remove gutta percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseal (obturate) with Hygenic gutta percha before preparing the post space.
4. To begin creating the post space, select the ParaPost Taper Lux Drill diameter which corresponds to the last Gates-Glidden drill used in Step 3. (See “Approximate Diameter Comparisons” chart, **Fig. A**). To help maintain the preplanned depth of the preparation, ParaPost Taper Lux Drills have been premarked to identify depths of 7, 9 and 11 mm from the apical end of the drill.

 ParaPost Taper Lux Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle. When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

5. Sequentially step up to the next larger ParaPost Taper Lux Drill until the preplanned diameter and depth are achieved **(Fig. 02)**.
6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box (Fig. 03).
7. Select the ParaPost Taper Lux Post that corresponds to the last drill used in the preparation of the post space. Remove colored band from post head. Insert post into post space **(Fig 04)**. Carefully check for occlusal clearance. Remove post; shorten as necessary from occlusal end with a Carborundum separating disk **(Fig. 05)**. Wear full face safety mask when shortening post. (Wet the post with water before cutting and rotate the post when cutting to ensure an even cut through each fiber.)
8. After try-in and before cementing, gently clean post with an alcohol wipe.
9. (Optional) Brush silane onto the post and allow to dry.
10. Condition the canal and/or apply primer/adhesive as per cement manufacturer's instructions (we recommend ParaCem Universal DC Cement).
11. Use the cement components according to manufacturer's instructions and apply to post surface. Use a Lentulo Spiral to coat the walls of the post space with cement, if setting time of cement will allow **(Fig. 06)**.
12. Seat post immediately, slowly inserting the post to full depth, allowing excess cement to vent **(Fig. 07)**. Apply pressure for about 60 seconds.
13. In the event light-curing is required, wipe away any excess cement before activating the light and follow manufacturer's instructions. Place curing light probe directly over the post, applying gentle pressure with the curing light tip. Light-cure for a minimum of 60 seconds (depending on volume and depth). If you are using an LED light, make sure your light-cured cement contains the Camphorquinone (CQ) Photoinitiator. If you have any questions, contact the cement manufacturer.
14. Complete the core and fabricate the final restoration **(Fig. 08)**.

Cleaning and Sterilization

Posts: Clean ParaPost Taper Lux post(s) in ultrasonic cleaner for 5 minutes. Remove and thoroughly rinse in tap water; dry with a paper towel. Place post(s) into a small, clean covered glass container. Add a ready-to-use FDA cleared liquid sterilant solution containing 7.50% hydrogen peroxide and 0.85% phosphoric acid (Sporox II**), covering the post(s). Let stand at room temperature for a minimum of 6 hours but not to exceed 24 hours. Pour off solution into sink and add sterile water to glass container covering the post(s). Gently agitate rinse solution for two minutes, then discard solution into sink and dump post(s) onto sterile sponge in clean container. Affix cover and store until ready for use.

Follow solution manufacturer's instructions.

Never re-use solution - never add fresh solution to used solution.

Always discard solution after each use!

ParaPost Taper Lux posts may also be steam sterilized (12 - 15 minutes at 130˚-135˚C). Follow manufacturer's recommendations. Posts must be sealed in a sterilization bag before use in autoclave.

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Steam sterilize at 132°C for 12 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization has been achieved, place autoclave bag, with the drills inside, in storage until ready for use.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by following instructions for drills above.

Warning Notice: DO NOT BEND POSTS!

Post Removal Technique

1. Take a radiograph of the tooth that requires post removal.
2. Carefully estimate length and diameter of ParaPost Taper Lux post to be removed. If possible, consult patient's record.
3. Remove any restorative or core material which may be in the way to gain direct straight-line access to the post.

4. Wear full face safety mask when removing post from tooth. Indent the center of the post using a # 1/2 round bur at high speed.
5. Carefully drill an initial channel into the center of the ParaPost Taper Lux Post using a Kodex twist drill (cat. # K95 or K97) at slow speed.
6. Carefully drill to full depth through the center of the post using an end cutting Tenax Starter Drill (cat # Te-DC-1, Ø.036"/0.90 mm) or Tenax Drill (cat # Te-A-11, Ø.043"/1.1 mm). Check radiographically at each stage.
7. Gradually work up to original diameter using ParaPost Taper Lux Drills.
8. Re-treat endodontically, if necessary.
9. Restore with new post/core and restoration.

Always verify steps via radiograph.

* References and/or literature available upon request.

**Sporox II is a registered trademark of Reckitt & Coleman, Inc.

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Made in France ©2006 P/N 86828C
Coltène/Whaledent Inc.
Coltène/Whaledent AG
Coltène/Whaledent GmbH + Co. KG