ParaBond®

DE	Gebrauchsinformation	2
EN	Instructions for Use	4
FR	Mode d'emploi	6
IT	Istruzioni per l'uso	8
ES	Instrucciones de uso	10
PT	Instrucões de Utilização	12
NL	Gebruiksaanwijzing	14
DA	Brugsanvisning	16
SV	Bruksanvisning	18
NO	Bruksanvisning	20
FI	Käyttöohje	22
EL	Οδηγίες χρήσης	24
TR	Kullanma Talimatları	26

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Instructions for Use

EN

Description
ParaBond Adhesive is a chemical cured, self-conditioning Adhesive System for enamel and dentin. It consists of a Non-Rinse Conditioner and a chemical-curing adhesive (Adhesive A and B).

Composition
ParaBond Non-Rinse Conditioner (NRC) contains:
Water
Acrylamidosulfonic acid
Methacrylate

ParaBond Adhesive A contains:
Methacrylates
Polyalkenoate
Initiators

ParaBond Adhesive B contains:
Ethanol
Water
Initiators

Areas of application
Adhesive bonding of chemical and dual curing cements and core build up materials.

Contraindications
Hypersensitivity due to any of the elements within ParaBond. Inadequate oral hygiene. If the working area can not be kept completely dry during application.

Side effects
Elements of ParaBond may lead to sensitivity or an allergic reaction for patients with a predisposed condition.

Interactions
Phenolics and other substances (e.g. zinc oxide eugenol) may not come in contact with ParaBond, since they will inhibit polymerization.

Application onto enamel and dentin
Conditioning the preparation/cavity
1. Clean the preparation/cavity with water and then dry excess water using a gentle stream of air for 2 s. Do not over dry the dentin.

Applying ParaBond Non-Rinse Conditioner
2. Dispense the Non-Rinse Conditioner into the mixing well.
3. Apply Non-Rinse Conditioner onto the entire preparation/cavity using a brush. Massage for 30 s.
4. Dry excess Non-Rinse Conditioner using a gentle stream of air for 2 s.

Alternative: 35 % phosphoric acid can be used for etching instead of ParaBond Non-Rinse Conditioner (according to the manufacturer's instructions for use).

Applying ParaBond Adhesive
5. Mix one drop of Adhesive A with one drop of Adhesive B.

Note: The working time is 2 min from the start of mixing (when material is stored in the refrigerator). Higher temperatures will accelerate the setting time of the material.
6. Apply mixed adhesive onto the preparation/cavity using a brush. Massage for 30 s.
7. Dry adhesive residue using a gentle stream of air for 2 s.

Important: The overall time between applying the Adhesive and the cementation should not exceed 5 min. If this time is exceeded, repeat the procedure again starting at 5. Excessive amounts of adhesive can accelerate the setting time of the cement.

Conditioning the inner surfaces/contact surfaces of the restoration
Always condition the contact surfaces of the restoration according to the manufacturer's instructions for use.

Note: After conditioning the restoration, rinse it thoroughly, dry with oil free air and protect against contamination (e.g. moisture, fingerprints). The bond strength can be improved by applying an extremely thin layer of adhesive to the inner surfaces of the restoration. This can however, accelerate the polymerization of chemical and dual cured cements.

Important: If the adhesive layer is too thick, it can impair the fitting accuracy of the final restoration. The adhesive can accelerate the setting time of the cement. Avoid application of excessive amounts.

Emergency measures
By direct contact with the oral mucosa, simply rinse using water. By contact with the eyes, thoroughly rinse with water (10 min) and consult an eye specialist showing the specialist these instructions for use.

Note
Close bottles after each use to prevent cross contamination. Keep the material out of reach of children! Supply only to dentists and dental laboratories or upon their instructions.

Caution
Federal law restricts this device to sale by or on the order of a dentist.

Shelf life and labelling
The expiration date and LOT number are printed on the containers and external packaging.
Do not use after the expiration date.

Storage
The system should be stored at a temperature of 4-8 °C / 39-46 °F.

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Coltène/Whaledent AG