FUSION DURALINK

PORTUGUÊS

Generalidades

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ENGLISH

Generalities

Fusion Duralink is an adhesive system of 4th generation, where primer and adhesive are in separate bottles, what makes the technique less susceptible to errors. Fusion Duralink is composite of a bottle of primer, a bottle of adhesive photopolymerizable and a bottle of chemical adhesive for cases where there is no possibility to reach the light.

Indications

Fusion Duralink is an adhesive system indicated for direct composite resin restorations and for the adhesive cementation of prefabricated posts (metal and fiber), cast metal posts and indirect restorations (metal, porcelain and composite resin). Its components can also be used for the bonding of orthodontic brackets and amalgam.

Composition

- Phosphoric acid: orthophosphoric acid; water; thickener; pigments
- Dentin primer: methacrylate esthers; ethanol; water
- Adhesive: methacrylate esthers; photoinitiators; polymerization accelerators
- Catalyst: methacrylate esthers; polymerization accelerators

Instructions for use

A. Procedure on ENAMEL only:

1. Certify that all conditions are adequate (isolation of operative field, perfect moisture control).
2. Apply the phosphoric acid gel on enamel for 30 seconds.
3. Rinse the acid gel with a water spray for 15 s.
4. Dry with air until the complete removal of moisture.
5. With a disposable brush, apply a thin and uniform layer of FUSION DURALINK ADHESIVO on the etched surface.

⚠ *Important: IT IS NOT NECESSARY TO DRY THE ADHESIVE.*

⚠ *In procedures on enamel only, it is not necessary to apply the primer.*

6. Light cure the adhesive for 10 s.
7. Continue the procedure.

B. Procedure on ENAMEL and DENTIN:

1. Certify that all conditions are adequate (isolation of operative field, perfect moisture control).
2. Apply the phosphoric acid gel on enamel for 30 s and on dentin for 15 s.
3. Rinse the acid gel with a water spray for 15 s.
4. Remove the excess of water with absorbent paper or a clean brush, leaving the surfaces slightly wet (shiny). If too dry, rewet them with a water spray and remove the excess of water again.

Important: DO NOT USE THE AIR JET TO REMOVE EXCESS OF WATER

5. With a disposable brush, apply a layer of FUSION DURALINK PRIMER on dentin and rub it slightly for 30 s.
6. Apply a light air jet from a distance of 10 cm for 10 s; the dentin should be shiny and without excess of the Primer; if there is a non-shiny spot, reapply the Primer and dry it again.
7. With another disposable brush, apply a thin and uniform layer of FUSION DURALINK ADHESIVE on enamel and dentin.

⚠ *Important: IT IS NOT NECESSARY TO DRY THE ADHESIVE*

8. Light cure the adhesive (ADHESIVE) for 10 s.
9. Continue the procedure.

C. Procedure with PORCELAIN (feldspathic or leucite reinforced):

1. To increase the adhesion to these porcelains, they must be etched with ANGELUS 10% HYDROFLUORIC ACID GEL for 1 minute.
2. After the gel is totally removed with a water spray, dry the porcelain with air and apply a layer of a silane (SILANO ANGELUS). Wait for 1 minute and apply a uniform layer of FUSION DURALINK ADHESIVE. Light cure for 10 s.

D. Cementation of intraradicular posts:

⚠ *Attention: no step during post cementation requires light curing.*

1. Certify that all conditions are adequate (isolation of operative field, perfect moisture control).
2. Apply the phosphoric acid gel in the root canal for 15 to 30 s.
3. Rinse the acid gel with a water spray for 15 s.
4. Remove the excess of water with absorbent paper points or a clean brush.
5. With a disposable brush, apply a layer of FUSION DURALINK PRIMER in the root canal and rub it slightly for20 s.
6. Apply a light air jet from a distance of 10 cm for 10 s; remove the excess of PRIMER with an absorbent paper point.
7. With another disposable brush, apply a thin and uniform layer of FUSION DURALINK SELF-CURE BOND (CATALYST) on the root canal.
8. Use resin cements that are dual cured or self cured only, such as ANGELUS CEMENT POST. Never use composite resins or resin cements that are light cured only.

⚠ *Important: Before cementation, the selected post must be prepared following these instructions:*

a. Clean it with alcohol.
b. Apply a layer of a silanating agent (SILANO ANGELUS), wait for 1 minute and dry with air.
c. Apply on the post a layer of FUSION DURALINK SELF-CURE (CATALYST).

Warnings

1. For transportation and storage of this product, temperatures should be kept between 21 e 24°C.
2. Do not use this product in patients with reported sensitivity to its components.
3. Use gloves, mask and eye protection when handling the product.
4. The curing light must have at least 400 mW/cm^2 of power.

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ATTENTION: This product must be used according to the instructions described in this manual. The manufacturer is not responsible for failure or damage caused by incorrect handling or use.

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