ParaCem® Universal DC

DE	Gebrauchsinformation	2
EN	Instructions for Use	8
FR	Mode d'emploi	14
ES	Información para el uso	20
IT	Informazioni per l'uso	26
SV	Bruksanvisning	32
NL	Gebruiksaanwijzing	38
DA	Brugsanvisning	44
NO	Bruksanvisning	50
FI	Käyttöohje	55
EL	Οδηγίες χρήσης	61
PT	Instrucões de Utilização	68
TR	Kullanma Talimatı	74

www.coltene.com

COLTENE®

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Instructions for use (English)

Description
ParaCem® Universal DC is a dual-curing adhesive system for cementing endodontic posts, core build-up, crowns, bridges, inlays, onlays and veneers. It consists of a Non-Rinse Conditioner, a chemically curing adhesive system (Adhesive A and Adhesive B) and a dual-curing cement (base and catalyst).

Composition
ParaBond® Non-Rinse Conditioner contains:
Water
Acrylamidosulfonic acid
Methacrylate

ParaBond® Adhesive A contains:
Methacrylates
Polyalkenoate
Initiators

ParaBond® Adhesive B contains:
Ethanol
Water
Initiators

ParaCem® Catalyst contains:
Methacrylates
Barium glass silanized
Amorphous silica
Initiators

ParaCem® Base contains:
Methacrylates
Barium glass silanized
Amorphous silica
Initiators

Areas of application
- permanent cementation of all types of endodontic posts, crowns, bridges, inlays and onlays (porcelain, metal and composite).
- core build-up
- Other applications: all-porcelain and composite veneers.

Contraindications
Allergies to the components of ParaCem® Universal DC. If the tooth cannot be kept completely dry during application. Patients with poor oral hygiene.

Side effects
There are no known systemic side effects. Contact allergies to products with a similar composition have been observed in isolated cases.

Interaction
Materials containing eugenol and clove oil can impair polymerization of ParaCem® Universal DC. Avoid using zinc oxide-eugenol products with ParaCem® Universal DC.

Application
General
As ParaCem® Universal DC is a dual-curing material, it is light sensitive and should not be exposed to blue light or strong lighting (surgical light, ambient light etc.) before polymerization. Protect pre-dispensed ParaCem® Universal DC, e.g. with an opaque covering (not blue), or dispense it immediately prior to application.

A. Posts
A1. Cementing Post followed by core build-up
1. Select a suitable root post (e.g. a ParaPost® post).
2. Prepare the root canal according to the instructions of the post manufacturer.

3. Select the appropriate coreformer and adjust with scissors.

Applying ParaBond® Non-Rinse Conditioner
4. Dispense the Non-Rinse Conditioner into the mixing well.
5. Apply the Non-Rinse Conditioner to the pretreated root canal and on the contact surfaces (preparation/cavity) for the core build-up with a brush and massage in for 30 s
6. Remove excess Non-Rinse Conditioner from the root canal with a paper point and then dry with a light jet of air.

Applying ParaBond® Adhesive
7. Mix one drop of Adhesive A with one drop of Adhesive B. Then apply the mixed Adhesive to the root canal and on the contact surfaces for the core build-up with a brush.
8. The working time is 2 min from the start of mixing (from the refrigerator). Higher temperatures accelerate setting. Leave the mixed Adhesive in the root canal and on contact surfaces for 30 s.
9. Remove excess Adhesive from the root canal with a paper point and then dry the entire bond layer with a light jet of air (Too much Adhesive residue accelerates the curing of the cement in the root canal).

Important: The time between applying the Adhesive and cementation should not exceed 5 min. If this time is exceeded, repeat the procedure from A.7.

Cementing the post
10. Dispense the base and catalyst in a ratio of 1:1 on the mixing pad and mix for 20–30 s until it forms a homogeneous paste.
11. The working time for the cement is 3 min from the start of mixing (from the refrigerator). Higher temperatures accelerate setting.
12. Coat the root post completely with the mixed cement and insert the post into the root canal using light pressure. If required, apply the cement into the canal using a lentulo. Start building up the core (e.g. ParaPost® ParaCore®) as soon as the cement has hardened (approx. 4 min after cementation).

The mixed cement can be light cured for 30 s to accelerate polymerization or to reduce the inhibition layer.

Application with ParaPost® ParaForm Coreformer
13. Dispense the base and catalyst in a ratio of 1:1 onto the mixing pad and mix for 20–30 s until it forms a homogeneous paste. The working time for the cement is 3 min from the start of mixing (from the refrigerator.) Higher temperatures will accelerate setting. *(Curing time is approx. 4 min).*
14. Fill the mixed Base and Catalyst directly into the coreformer. Apply the coreformer to the tooth and position with slight pressure. Remove any remaining excess with a suitable instrument (see removing excess material).

To accelerate polymerization, ParaCem® can be light cured (e.g. Coltène® Coltolux) for 40-60s (20s per side/surface at 2mm layer thickness).

15. Remove the coreformer using a scalpel or rotary instrument. Prepare the core reconstructed with ParaCem® in the usual way.

A.2 Post Only
Proceed as described A.1 section 1-12 without select coreformer and application of ParaBond® on the contact surfaces for core build-up.

B. Core Build-up
Pre-treatment according to the instructions of the manufacturer.

Conditioning the preparation/cavity
1. Clean the preparation/cavity with water and then dry the excess

water with a light jet of air (2 s). Do not over-dry the dentin.

Applying ParaBond® Non-Rinse Conditioner

2. Dispense the Non-Rinse Conditioner into the mixing well.
3. Then apply the Non-Rinse Conditioner to the entire preparation/cavity with a brush and massage in for 30 s
4. Dry the excess Non-Rinse Conditioner with a light jet of air.

Alternative: 35% phosphoric acid can be used for etching instead of ParaBond® Non-Rinse Conditioner.

Applying ParaBond® Adhesive

5. Mix one drop of Adhesive A with one drop of Adhesive B.
6. The working time is 2 min from the start of mixing (from the refrigerator). Higher temperatures will accelerate setting.
7. Then apply the mixed Adhesive to the preparation/cavity with a brush and massage in for 30 s.
8. Dry the excess Adhesive with a light jet of air.

Important: The time between applying the Adhesive and cementation should not exceed 5 min If this time is exceeded, repeat the procedure from B.5. Avoid pools of Adhesive, as these can impair the fit of the restoration.

Application with ParaPost® ParaForm Coreformer

9. Select the appropriate coreformer and cut to size with scissors.
10. Dispense the base and catalyst in a ratio of 1:1 onto the mixing pad and mix for 20–30 s until it forms a homogeneous paste. The working time for the cement is 3 min from the start of mixing (from the refrigerator.) Higher temperatures will accelerate setting. *(Curing time is approx. 4 min).*
11. Fill the mixed Base and Catalyst directly into the coreformer. Apply the coreformer to the tooth and position with slight pressure. Remove any remaining excess with a suitable instrument (see removing excess material).

To accelerate polymerization, ParaCem® can be light cured (e.g. Coltène® Coltolux) for 40–60s (20s per side/surface at 2mm layer thickness).

12. Remove the coreformer using a scalpel or rotary instrument. Prepare the core reconstructed with ParaCem® in the usual way.

C. Crowns/Bridges, Inlays/Onlays

Conditioning the preparation/cavity

1. Clean the preparation/cavity with water and then dry the excess water with a light jet of air (2 s). Do not overdry the dentin.

Applying ParaBond® Non-Rinse Conditioner

2. Dispense the Non-Rinse Conditioner into the mixing well.
3. Then apply the Non-Rinse Conditioner to the entire preparation/cavity with a brush and massage in for 30 s
4. Dry the excess Non-Rinse Conditioner with a light jet of air.

Alternative: 35 % phosphoric acid can be used for etching instead of ParaBond® Non-Rinse Conditioner.

Applying ParaBond® Adhesive

5. Mix one drop of Adhesive A with one drop of Adhesive B.
6. The working time is 2 min from the start of mixing (from the refrigerator). Higher temperatures accelerate setting.
7. Then apply the mixed Adhesive to the preparation/cavity with a brush and massage in for 30 s.
8. Dry the excess Adhesive with a light jet of air.

Important: The time between applying the Adhesive and cementation should not exceed 5 min If this time is exceeded, repeat the procedure from C.5. Avoid pools of Adhesive, as these can impair the fit of the restoration.

Conditioning the inner surfaces/contact surfaces of the restoration

Always condition the contact surfaces of the restoration according to the manufacturer's instructions.

Note: After conditioning the restoration, rinse it thoroughly, dry it with oil free air and protect against contamination (e.g. moisture, fingerprints). The bond strength can be improved by applying an extremely thin layer of Adhesive A/B to the inner surfaces of the restoration. This can, however, accelerate polymerization of the ParaCem® Universal DC.

Important: If the Adhesive layer is too thick, it can impair the fit of the restoration.

Cementing crowns / bridges / inlays / onlays

1. Dispense the base and catalyst in a ratio of 1:1 onto the mixing pad and mix for 20–30 s until it forms a homogeneous paste.
2. The working time for the cement is 3 min from the start of mixing (from the refrigerator.) Higher temperatures will accelerate setting. *(Curing time is approx. 4 min after cementation.)*
3. Apply cement with a suitable instrument (e.g. brush, spatula) to the inner surfaces of the restoration and/or if necessary (with concave shapes to prevent trapped air bubbles) to the preparation.
4. Then place the restoration in position using slight pressure.
5. Remove rough excess (e.g. with a brush, spatula) and then hold the restoration in position with increased pressure.

The cement can be light cured for 30 s to accelerate polymerization or to reduce the inhibition layer. (Light-curing is not necessary for the material to harden completely).

D. Veneers

Conditioning the preparation

1. Clean the preparation with water and then dry the excess water with a light jet of air (2 s). Do not over-dry the dentin.

Applying ParaBond® Non-Rinse Conditioner

2. Dispense the Non-Rinse Conditioner into the mixing well.
3. Then apply the Non-Rinse Conditioner to the entire preparation with a brush and massage in for 30 s
4. Dry the excess Non-Rinse Conditioner with a light jet of air.

Alternative: 35% phosphoric acid can be used for etching instead of ParaBond® Non-Rinse Conditioner.

Applying ParaBond® Adhesive

5. Mix one drop of Adhesive A with one drop of Adhesive B.
6. The working time is 2 min from the start of mixing (from the refrigerator). Higher temperatures will accelerate setting.
7. Then apply the mixed Adhesive to the preparation with a brush and massage in for 30 s.
8. Dry the excess Adhesive with a light jet of air.

Important: The time between applying the Adhesive and cementation should not exceed 5 min If this time is exceeded, repeat the procedure from D.5. Avoid pools of Adhesive, as these can impair the fit of the restoration.

Conditioning the inner surfaces/contact surfaces of the restoration

Always condition the contact surfaces of the restoration according to the manufacturer's instructions.

Note: After conditioning the restoration, rinse it thoroughly, dry it with oil free air and protect against contamination (e.g. moisture, fingerprints). The bond strength can be improved by applying an extremely thin layer of Adhesive A/B to the inner surfaces of the restoration. This can, however, accelerate polymerization of the ParaCem® Universal DC.

Important: If the Adhesive layer is too thick, it can impair the fit of the restoration.

Cementing Veneers
Recommendation
1. Dispense the base and catalyst in a ratio of 1:1 onto the mixing pad and mix for 20–30 s until it forms a homogeneous paste.
2. The working time for the cement is 3 min from the start of mixing (from the refrigerator.) Higher temperatures accelerate setting. *(Curing time is approx. 4 min after cementation.)*
3. Apply cement with a suitable instrument (e.g. brush, spatula) to the inner surfaces of the restoration and/or if necessary (with concave shapes to prevent trapped air bubbles) to the preparation.
4. Then place the restoration in position using slight pressure.
5. Remove rough excess (e.g. with a brush, spatula) and then hold the restoration in position with increased pressure.

The cement can be light cured for 30 s to accelerate polymerization or to reduce the inhibition layer. *(Light-curing is not necessary for the material to harden completely).*

Optional: ParaCem® Bases A2 and B3 can be used as light-curing cement, If there is sufficient light transmission through the complete restoration (according to the instruction of the manufacturer).
We do not recommend polymerizing ParaCem® base white by light-curing only because of the increase in opacity. A mixture of base and catalyst should always be used in this case to ensure full polymerization.

Removing excess material
1. While holding the restoration in position, remove any remaining excess with a fine spatula, brush or other suitable instrument. (The use of matrices, wedges, rubber dams etc. shall prevent ParaCem® Universal DC from flowing into the sulcus and interproximal region, as it is very difficult to remove after it hardens.)
2. As soon as the consistency of ParaCem® Universal DC becomes gel-like, do not remove any more excess. The inhibition layer (contact of the outer layer of cement with the oxygen in the air) that always forms during chemical polymerization at the margins of the restoration can be compensated for with a very thin amount of excess material.

Preparing and polishing
Hardened excess can be removed with a scaler, very fine diamonds or flexible discs. Finishing and polishing strips can be used in the interproximal region. Use silicone rubber polishers and brushes for finishing.

Emergency measures
On direct contact with the oral mucosa simply rinse with water. On contact with the eyes, thoroughly rinse with water (10 min) and then consult an eye specialist and show the specialist these instructions for use.

Note
Close syringes and bottles after each use to prevent any mutual contamination.
Keep the material out of reach of children!
Supply only to dentists and dental laboratories or upon their instructions.

Caution
Federal law restricts this device to sale by or on the order of a dentist.

Shelf life and labelling
The expiration date and LOT number are printed on the containers and external packaging. Do not use after the expiration date.

Storage
The system should be stored at a temperature of 4–8 °C / 39–46 °F.

Date of issue
05-2012

For MSDS see
www.coltenewhaledent.com

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Coltène/Whaledent AG