RealSeal Self-Etch™ System

Index

INSTRUCTIONS FOR USE

Hrvatski/Croatianp2
Česky/Czechp4
Dansk/Danishp6
Nederlands/Dutchp8
Englishp10
Eesti/Estonianp12
Suomi/Finnishp14
Français/Frenchp16
Deutsch/Germanp18
Elinika/Greekp20
Magyar/Hungarianp22
Íslenska/Icelandicp24
Italiano/Italianp26
Lietuviška/Lithuanianp28
Norsk/Norwegianp30
Polski/Polishp32
Português/Portuguesep34
Română/Romanianp36
Srpski/Serbianp38
Slovenčina/Slovakianp40
Slovenščina/Slovenianp42
Español/Spanishp44
Svenska/Swedishp46
Türkçe/Turkishp48

SybronEndo

077-3134 Rev. A

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English

INSTRUCTIONS FOR USE:

The *RealSeal* Obturation System has been proven to result in obturations that resist leakage significantly better than gutta percha-based obturations with the added benefit of strengthening the root. The centerpiece of this system is a soft resin obturation material called Resilon®*. *Resilon* material is a formulation of polymers of polyester with unique fillers and radiopacifiers in a soft resin matrix. *Resilon* material offers all the features gutta percha provides with the added benefits of enhanced sealing to minimize microleakage and the ability to thermoplasticize at lower temperatures.

The *RealSeal Self-Etch* System with *Resilon* obturation material:

- Is available in both points and pellets
- Is highly radiopaque
- Looks and handles like gutta percha
- Is fully biocompatible and easily retreatable

RealSeal Self-Etch dual-cure, hydrophilic resin sealer bonds to both *RealSeal* obturators and to dentin in the canal to create a monoblock. An immediate coronal seal is formed by light curing the coronal aspect. With obturation material in both points and pellets, the *RealSeal* System works effectively and seamlessly with all endodontic filling techniques.

Important Notes:

- Complete cleaning and shaping of the canal is essential.
- If calcium hydroxide has been used as an intracanal medicament, completely remove it with EDTA before obturation.
- Always apply additional sealer after downpacking when using the vertical condensation technique to ensure that the sealer completely lines the walls of the canal.
- Use alcohol, instead of chloroform, to remove excess sealer from the chamber floor.
- After light curing the coronal aspect for an immediate coronal seal, place a bonded composite restoration to complete the treatment. *RealSeal Self-Etch* Sealer will self cure within 45 minutes to seal the entire canal.

Indications:

- The *RealSeal* System is indicated for the obturation of root canals.

Contraindications:

- *RealSeal* points and pellets are contraindicated for use with sealers formulated for use with gutta percha.

Storage and Use:

1. *RealSeal* Points and Pellets – Store and use at room temperature, 65˚-75˚F (18˚-24˚C).

2. *RealSeal Self-Etch* Sealer – Refrigeration required. Store at 35˚-55˚F (2˚–12˚C). Use at room temperature.

Instructions for Use

1. Irrigate the prepared canal using the following sequence:

a. Cleanse the canal with sodium hypochlorite (NaOCl) and dry.

b. Remove the smear layer with EDTA.

Note: If EDTA with a surfactant is used, flush the canal with sodium hypochlorite to remove residual surfactant which can interfere with bonding.

c. Rinse the canal with sterile water or 2% chlorhexidine and dry.

Important Irrigation Notes: Sodium hypochlorite and peroxide-based lubricants may affect bonding, but are essential for cleaning during the irrigation process. To ensure the best possible bond, please observe the following recommendations:

a. Sodium hypochlorite (NaOCl) is a required irrigant, but must not be the final irrigant used in the canal.

b. Dry the canal with paper points. Do not dessicate the canal with alcohol.

2. Select the appropriate sizes of *RealSeal* points for vertical or lateral condensation techniques and/or pellets for backfilling.

3. Sterilize the obturation materials. Place *RealSeal* points and/or pellets into 2% chlorhexidine or other suitable disinfectant solution for 60 seconds. Rinse with sterile water and dry before use.

Note: Do not use sodium hypochlorite to disinfect *Resilon* points or pellets.

4. Trial seat *RealSeal* points to ensure the selected point(s) reach the working length with tugback. Take a radiograph to confirm apical positioning.

5. Dispense *RealSeal Self-Etch* Sealer onto a mixing pad.

6. Apply sealer to the walls of the canal using any of the following techniques:

i. A master cone coated with sealer

ii. A paper point coated with sealer

(Continued on the next page.)

English

INSTRUCTIONS FOR USE:

iii. A small hand file with length control, rotated counter clockwise within the canal

iv. A lentulo spiral, 3mm from the apex at 300 RPM

Note: Sealer may also be dispensed directly into the canal by using the intra-oral tips supplied.

7. Fill the canal with *RealSeal* points or pellets using your normal, preferred technique.
Note: If using the vertical condensation technique, apply additional sealer after downpacking. The intra-oral tips supplied with *RealSeal Self-Etch* Sealer are ideal for this use. If using lateral condensation, coat each accessory point with sealer.

8. Remove any excess sealer from the chamber floor with alcohol. Do not use chloroform.

9. Take a check radiograph to ensure that the canal is completely filled.

10. Light cure the coronal surface of the obturation for 40 seconds to create an immediate coronal seal. *RealSeal Self-Etch* Self-etch Sealer will self-cure in the canal within 45 minutes creating a monoblock that effectively resists leakage.

11. Seal the access cavity.

Note: Although *RealSeal* obturations create an immediate coronal seal when light cured and seal the canal within 45 minutes upon self cure, the placement of a bonded composite over the pulpal floor is recommended. The patient should be encouraged to receive the final restoration as soon as possible. The quality of the restoration is a significant factor in the success of the endodontic treatment.

Backfilling with Heated Delivery Systems

If using Obtura™* or System B™** to backfill the obturation with *RealSeal* Hi Flow Pellets, use the following settings:

• Obtura™†:

- **o** For the 25 gauge needle tips, set the temperature to 125°C.
- **o** For the 23 gauge needle tips, set the temperature to 115°C.
- **o** For the 20 gauge needle tips, set the temperature to 105°C.

• System B: Set the temperature to 150°C and the power to 10.

† These temperature settings are initial recommendations when using Obtura™* and may require slight adjustment depending upon the calibration of the unit.

* *Obtura* is a trademark of Obtura/Spartan.

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Resilon is a registerd trademark of Pentron Clinical Technologies, LLC.
** System B is a trademark of Sybron.