


GUTTA - CONDENSOR

Ref A0242 - A 0244

P **Português**	Este manual de utilização está também disponível, quando requisitado, nas seguintes línguas: português, holandês, dinamarquês, sueco, finlandês e grego.
NL **Nederlands**	Deze gebruiksaanwijzing is, op aanvraag, eveneens verkrijgbaar in devolgende talen: Portugees, Nederlands, Deens, Zweeds, Fins en Grieks.
FIN **Suomi**	Käyttöohje on saatavana myös seuraavilla kielillä: portugali, hollanti, tanska, ruotsi, suomi ja kreikka.
S **Svenska**	Denna bruksanvisning finns även att tillgå på följande språk:portugisiska, tyska, danska, svenska, finska och grekiska.
DK **Dansk**	Denne brugsanvisning kan også rekvireres på følgende sprog: portugisisk, hollandsk, dansk, svensk, finsk og græsk.
GR **Ελληνική.**	Αυτή η οδηγία είναι επίσης διαθέσιμη, κατόπιν ζήτησης, στις ακόλουθες γλώσσες: Πορτογαλέζικη, Γερμανική, Δανέζικη, Σουηδική, Φιλανδέζικη και Ελληνική.

Visit our website: www.dentsplymaillefer.com

FOR DENTAL USE ONLY **GB**

DIRECTIONS FOR USE

1) INDICATIONS FOR USE

These instruments are to be used only in a clinical or hospital environment, by qualified users.

Application field: The Gutta-Condensor is used to mechanically condense gutta percha in the root canal, after instrumentation, cleaning and disinfection are finished.

2) CONTRAINDICATIONS

- The Gutta-Condensor must not be used in curved canals, as in this clinical situation the risk of separation of the instrument is high.
- The Gutta-Condensor is not recommended for use in patients who have a known allergy to nickel. Use of this product in these patients may cause: difficulty breathing, swollen face or eyes, hives or rash. Should any of these symptoms occur, the patient should be advised to contact his/her dental professional immediately.

3) WARNINGS

- The use of the Gutta-Condensor is not risk-free. The main danger consists in a possible penetration of the gutta percha into the periapical region. The following reasons may be responsible for such a case: the use of a too-small or too-fine master cone, a too-small Gutta-Condensor, excessive speed or axial pressure, or too-deep insertion of the Gutta-Condensor in the root canal. Even a small amount of gutta percha in the periapical region may cause severe clinical reactions. There is no resorption to be expected from the periapical macrophages.
- For multicuspid mandibular teeth, an overfilling reaching the inferior dental nerve may provoke irreversible damage, such as diasthesia, paraesthesia, anaesthesia, or paralysis with a possible labiodynisis.
- The Gutta-Condensor is not recommended for use in patients who have a known allergy to nickel. Use of this product in these patients may cause: difficulty breathing, swollen face or eyes, hives or rash. Should any of these symptoms occur, the patient should be advised to contact his/her dental professional immediately.

4) PRECAUTIONS

- Never use the Gutta-Condensor in counter-clockwise rotation (danger of getting the stuck in the canal or in the bone).
- Do not penetrate the tip of the instrument into the canal to a distance of less than 2 mm from the apex.
- Do not oppose strong resistance when the Gutta-Condensor is backing out of the canal

5) ADVERSE REACTIONS

- For multicuspid mandibular teeth, an overfilling reaching the inferior dental nerve may provoke irreversible damage, such as diasthesia, paraesthesia, anaesthesia, or paralysis with a possible labiodynisis.
- The Gutta-Condensor is not recommended for use in patients who have a known allergy to nickel. Use of this product in these patients may cause: difficulty breathing, swollen face or eyes, hives or rash. Should any of these symptoms occur, the patient should be advised to contact his/her dental professional immediately.

6) STEP-BY-STEP INSTRUCTIONS

See scheme above

1. Establish the working length and apical stop of the root canal.
2. Insert the gutta percha point and the Gutta-Condensor into the canal, and press the Gutta-Condensor against the gutta percha.
3. Descend the Gutta-Condensor into the canal, rotating in a clock-wise direction at 8'000 min-1. Do not penetrate the Gutta-Condensor to a distance of less than 2 mm from the apex.
4. Back the Gutta-Condensor out of the canal, without opposing strong resistance.

7) DISINFECTION, CLEANING AND STERILIZATION

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GB **Reprocessing procedure for dental instruments and implantable radicular devices.**

Foreword

For hygiene and sanitary safety purposes, all instruments not marked "sterile" must be cleaned, disinfected and sterilized before each usage to prevent any contamination according to ISO 17664. This concerns the first use as well as the subsequent ones.

Area of application

A1. Instruments:
Cutting instruments, (hand and engine driven) such as:
- Endodontics instruments (files, broaches, reamers, enlargers, endodontic burs);
- Rotative cutting instruments (Diamond burs, tungsten carbide burs, stainless steel drills, carbon steel burs, reamers).
- Root canal filling instruments (Plugger, spreader, compactor);
- Supports, kits and systems of instrumental organization;
- Hand instruments.

A2. Implantable devices: Dentinal and radicular posts made of steel, titanium and fiber posts. Supports, kits and organization systems for the posts.

B. Filling material: Gutta percha, Thermafil obturation devices.

Exclusion

Paper points: Sterile paper points are recommended
Equipments as Motor, Contra-angle, Apex locator, where reprocessing procedures are included in the individual Direction for Use.
MTA, Glyde, Topseal.

General recommendation
- Single use marked instruments are not approved for re-use.
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments.
- For your own safety, please wear personal protective equipments (gloves, glasses).
- Use only disinfecting solution which is approved for their efficacy (DGHM-listing, CE marking, FDA approval)
- Tungsten Carbide burs, Plastic Support, Hand instruments and NiTi instruments are degraded by Peroxide hydrogen (H_2O_2) solution.
- NiTi Instruments are degraded if immerged more than 5 minutes in a solution of NaOCl at more than 5%.
- Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt. Do not use acid (pH < 6) or alkaline (pH > 8) solutions.
- After 5 cycles, degradation of marking may appear.

Step-by-step procedure

A. Devices

	Operation	Operating mode	Warning	A1. Instruments – First usage	A1. Instruments – Following usages	A2. Implantable devices
1.	Pre-Disinfection or Decontamination	- Soak immediately just after usage all instruments in a disinfectant solution combined with proteolytic enzyme if possible.	- Follow instructions and observe concentrations and immersion time given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments). - The disinfectant solution should be aldehyde free (to avoid blood impurities fixation). - Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (See general recommendation). - For visible impurities are observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.		x	
2a.	Automated Cleaning Disinfection	- Disassemble the devices (silicone stops to be removed) - Place them in a kit, support or container - Put them in the thermodisinfector (for at least 10 min at 93°C or Ao value > 3000).	- Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the thermodisinfector. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendation). - The thermodisinfector is not recommended for instruments made of aluminium, tungsten carbide or carbon steel. In case of chemical disinfection, there is a danger of remnants of the disinfectant on the instruments. - Sufficient rinsing step should be available in the program. Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata. - Use filtered air for drying the instruments or posts. - Use only approved thermodisinfector according to EN ISO 15883, maintain and calibrate it regularly. - If possible, an automated procedure is preferred.	x	x	x
2b.	Manual Cleaning Disinfection	- Disassemble the devices (silicone stops to be removed) - Immerse them in the disinfectant solution assisted with ultrasonic device if suitable. - Rinse the devices thoroughly with clean, demineralised or distilled water and then dry them with filtered compressed air.	- No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Avoid any contact between instruments or posts when placing in the solution, use kits, supports or container. - Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendation). - Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata. - If a disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before autoclaving.	x	x	x
3.	Inspection	- Inspect devices and sort out those with defects. - Assemble the devices (stops)	- Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument or posts. - Protect carbon steel bur with corrosion inhibitor before packaging.		x	x
4.	Packaging	- Pack the devices in "Sterilisation pouches".	- Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.	x	x	x
5.	Sterilization	- Steam sterilisation at: 134°C (237°F) during 18 min for endodontics instruments and root canal filling instruments 134 °C (237°F) during at least 5 min for all other devices.	- The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use fractionated vacuum or gravity (less preferred) autoclaves (according to EN 13060, EN 285). - Use validated sterilisation procedure according to ISO 17665-1 - Respect maintenance procedure of the autoclave device given by the manufacturer. - Use only the listed sterilization procedures.	x	x	x
6.	Storage	Keep devices in sterilization packaging in a dry and clean environment	- Sterility cannot be guaranteed if packaging is open, damaged or wet (check the packaging before using the instruments).	x	x	x

B. Filling material

1.	Disinfection	- Immerse the obturation devices in alcohol during 30 sec at ambient temperature.	- Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (See general recommendation).

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