


ADVANCED
TISSUE-MANAGEMENT
The easy way to aesthetics.
PARASORB Sombrero®
EN
RESORBA®
REPAIR AND REGENERATE

PARASORB Sombrero® – UNIQUE

The combination of two tried and tested products from RESORBA®: Membrane (RESODONT®) and Collagen Cone.


The Membrane:


- the barrier function of RESODONT® prevents ingrowth of connective tissue, whilst the surface structure of the membrane still permits growth of bone-forming cells.
- allows rapid epithelisation above the membrane
- saliva-proof closure protects against infections



The Collagen Cone:

- stabilisation of the blood clot
- stabilisation of the alveolus and thus maintenance of the sensitive vestibular lamella
- volumetric stability of the alveolar extension
- osteo- and angioconductive properties
- stable, bony healing of the alveolar socket
- improved long-term aesthetic and functional results






Experience with use of PARASORB Sombrero®


»Experience to date with the use of Sombrero has been very good. Both the adaptation to individual situations and the handling of the product on the patient are without problems.


Its use leads to reliable bony healing of the alveolus with a stable epithelial cover as the prerequisite for a subsequent implantation, maxillofacial tooth movement or prosthetic care.«

Prof. J. Kleinheinz,
Clinic and Policlinic for Oral and Maxillofacial Surgery
Muenster University Hospital, Germany

CONVINCING ARGUMENTS REGARDING USE







The combination:

- simple handling, as the membrane and cone are firmly attached to each other
- offers a defined regeneration matrix
- formation of a stable bony implant bed after a relatively short time
- only minimal detachment of the marginal gingiva requested at the alveolar ridge (minimises trauma / tissue supply largely maintained)
- improved long-term aesthetic and functional results
- unproblematic healing phase
- leads to a stable epithelial coverage
- completely absorbable

The indications:

- to cover and simultaneously fill the alveolus (socket preservation)
- to cover and simultaneously fill defects in guided bone regeneration (GBR)
- treatment of the alveolus and other bone defects

Sizes

PARASORB Sombrero®, sterile, is available in the following pack sizes:

REF	PU
RD0001	1 Membrane-Cone



RESORBA® REPAIR AND REGENERATE